

# Arvo™

Family of Receiver-in-Canal Products





## **CHOOSE ARRO™ FOR COMFORT, CLARITY AND TO ENRICH YOUR QUALITY OF LIFE.**

---

Hearing loss can diminish the quality of your life and create awkward moments in family, social and business settings. The Arro™— Receiver in Canal hearing instrument eliminates these uneasy moments. Placing the receiver in the ear canal opens your world to a natural listening experience without the clutter of background noise.

The advanced technology of the Arro™ provides clear hearing in most listening environments allowing you to enjoy dining out, talking with friends and listening to the soft sounds of nature. With the receiver in the ear canal, the “echo chamber” feeling common to other types of hearing aids is virtually eliminated. The open fitting allows your ear canal to remain clear and provides digital clarity without the distracting plugged-up feeling.

Arro™ miniature design gives you a discreet and personal solution for better hearing.

## **P2I NANO-COATING PROTECTS HEARING INSTRUMENTS**

---

Moisture, body oils, and ear wax can all contribute to hearing aid breakdown and failure. Audina offers protection from these elements using P2i nano-coating. This durable ultrathin fluoropolymer coating is applied to all internal and external components, acting as a water and oil repellent barrier. This hi-tech coating system ensures that the sensitive electronic components are protected to ensure optimal durability and performance. Ask your hearing specialist about P2i protection.

## ARRO™ MAJOR FEATURES


---

- Provides customized, natural comfort for a wide variety of hearing losses.
- Three (3) Hi-Def Receiver unit choices – Wideband, Power and High Power, allows your hearing professional to select the appropriate receiver for your specific type of hearing loss.
- 312 Battery for up to 180 hours of battery life.
- Adaptive Directionality automatically adjusts to noise in the environment, switching between microphone settings to adapt to the source of speech. This feature is available on select Arro™ models.
- Multi Memory Push Button – Access to multiple listening environments (up to 4) are available, such as restaurant dining, business conversations and musical concerts.
- Revolutionary Adaptive Feedback Cancellation System with constant monitoring - virtually eliminates whistling, known as feedback.
- Arro's advanced circuitry reduces annoyance and fatigue of background noise while giving clarity and low distortion to speech input, resulting in full and clear speech production.
- Environmental Recognition System (ERS) offers advanced reduction of environmental noise improving listening in the most extreme situations.
- Programmable Power-on Delay gives you time to properly place your Arro™ instrument without the normal feedback common when positioning a hearing instrument on or in your ear.

## THE IMPORTANT ADVANTAGES OF BINAURAL (TWO-EAR) FITTING

---

- A more natural “stereo sound” and an overall balanced feeling
- Improved ability to identify the direction of the sound, which is crucial for safety
- Hearing is less tiring, and listening is more pleasant
- Increased ability to understand speech in normal conditions and noisy situations
- Less volume needed, which reduces the interference of background noise and the chance of feedback
- Keeps both ears active, reducing the chance of the unaided ear losing its ability to hear or understand
- Greater comfort when loud noises occur
- Smoother tone quality and improved clarity



**THERE IS AN ARRO  
TO MEET EVERY LOSS**



For over twenty five years, our professional associates at Audina® Hearing Instruments, Inc. have designed and produced state-of-the-art hearing instruments to improve the lives of hearing loss patients.

- Privately owned American company and industry leader
- Comprehensive product line to accommodate all hearing loss needs
- Internationally recognized – exporting to over 45 countries

